

<b>Brand Name of Product</b>	Saliva Test Soil
<b>Generic Name of Product</b>	Saliva Test Soil, Modified Saliva Test Soil
<b>Product Code Number(s)</b>	STS-100ML, STS-NHC-45ML
<b>Intended Use</b>	Simulated use test soil for cleaning validations of oral devices.
<b>Range of Applications for Product</b>	For simulated use soiling of items for cleaning validations for oral devices.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>STS-100ML- Mixture of purified bovine proteins (hemoglobin, albumin), artificial saliva, mucin.</li> <li>STS-NHC-45ML- Mixture of purified bovine albumin, artificial saliva, mucin. No hemoglobin</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	<ul style="list-style-type: none"> <li>Store the separate vials with dry and wet components at room temperature.</li> <li>Keep away from light and heat.</li> </ul>
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	3 months from the date of manufacture

Instructions for Using Product	
<b>Description of Use(s)</b>	N/A
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	N/A
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>Add the provided solution S to the dry soil.</li> <li>Dissolve the components by shaking for 5 minutes. A laboratory vortex mixer may be used. Use within approximately 15 minutes of hydrating.</li> <li>Contaminate the item in a manner that reflects real-life soiling from use.</li> </ol>
<b>Interpretation of Test Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li><b>Hazard Information</b> <ul style="list-style-type: none"> <li>May be harmful if ingested or inhaled.</li> <li>May cause skin irritation.</li> <li>May cause eye irritation.</li> <li>May be irritant to mucous membranes and upper respiratory tract.</li> </ul> </li> <li><b>First Aid Measures</b> <ul style="list-style-type: none"> <li>After inhalation: Fresh air.</li> <li>After skin contact: Wash off with plenty of water. Remove contaminated clothing.</li> <li>After eye contact: Rinse with plenty of water with the eyelid held wide open. Call an Ophthalmologist.</li> <li>After swallowing: Make victim drink plenty of water. Call a physician.</li> </ul> </li> <li><b>Accidental Release Measures</b> <ul style="list-style-type: none"> <li>After spillage: Dilute spilled substance with plenty of water and absorb with absorbent material.</li> </ul> </li> </ul>
<b>Disposal</b>	Dispose of in accordance with regional and/or national regulations.

Reprocessing Instructions	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A

<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	N/A
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	N/A
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark, A Getinge Company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://www.artificialtestsoil.com">www.artificialtestsoil.com</a>